


# One Of The Smallest, But...THE STRONGEST!

## Introducing The Cure Cordless By Spring Health.

### Key Features of The Cure Cordless (TC-CL)

- Dual Switches: Ease of operation in any position.
- Digital timer with 8 settings form 2 to 20 seconds.
- High / Low output control from 600 to 1100 mw/cm<sup>2</sup> for step curing.
- Includes: Charging base stand with power supply, 8 mm lens, lens wrench, 2 light shields and 50 protective sleeves.



**Six Second Cure:** The CURE is available in a cordless model (TC-CL) and an ultra light (50 gram) model with 9 volt power supply and cord (TC-01). Both are highly rated by the leading evaluators of dental products. With an excellent rating in light intensity and depth of cure tests from Reality Publishing, Dental Advisor and a leading independent clinical research firm, the CURE offers fast, reliable curing in as little as 6 seconds on most modern composites (A3 or lighter shades).

**Output Never Degrades:** The CURE Cordless battery provides consistent output throughout the entire range of the charge.

**Light Weight/Small Size:** The CURE is manufactured with a high-tech, light weight, all metal housing making it more durable than units with plastic housings. Its small size resembles a high speed handpiece making it one of the smallest on the market yet one of the strongest.

**The Spring Health Guarantee:** Both models come with a one year warranty. Spring Health provides an additional satisfaction guarantee for 30 days. If you're not satisfied with the results you get from the CURE, simply return it to your dealer for a full refund.

Manufactured By:



## The Cure Original.

### Key features of The Cure (TC-01) (corded model)

- Dual Switches: Ease of operation in any position.
- High / Low output from 600 to 1100 mw/cm<sup>2</sup>.
- 16 second timed cycle, beeps every 2 seconds.
- Included with the TC-01 ORIGINAL: AC adapter power supply, 8 mm lens, lens wrench, 2 light shields, 50 protective sleeves, handpiece hanger.

*Also manufactured with a high-tech, light weight, all metal housing making it more durable than units sold with plastic housing.*



*Comes with wall-mounted handpiece hanger.*

