

Single-use FG Diamonds (10 & 25 Packs Only)



	Occlusal Reduction & Finishing		Wheel	Needle	Mosquito	Bulk Reducer	Depth Cutter	Flame	
Shape No.	368	811	909	889	392	6055	834	879	862
Size (1/10 mm)	023	033	035	008	016	018	021	018	012
Length (mm)	5.0	4.2	1.1	3.0	5.0	8.0	6.0	9.0	8.0
Spring No.	285.5	230	860	201.3	392-016	ST8	834-021	256.9	260.8
Supercoarse	000-0000	000-0000	000-0000	-	-	000-0000	-	000-0000	-
Coarse	-	-	-	-	-	-	-	-	-
Medium	-	-	-	000-0000	000-0000	-	000-0000	-	000-0000
Fine	000-0000	-	-	-	-	-	-	000-0000	000-0000

	Flame			Beveled Cylinder	Round			Inverted Cone			Double Inverted Cone	Flat End Cylinder	
Shape No.	863	863	863	878	801	801	801	805	805	805	813	835	836
Size (1/10 mm)	012	014	018	012	010	014	023	010	014	016	014	008	012
Length (mm)	10.0	8.0	9.0	8.0	1.0	1.4	2.3	1.0	2.0	2.0	2.0	4.0	5.0
Spring No.	260.10	265.8	270.9	250.8	115	120	135	312	313	320.2	1220	514.4	515.5
Supercoarse	-	-	-	-	-	-	000-0000	-	-	-	-	-	000-0000
Coarse	000-0000	000-0000	000-0000	000-0000	000-0000	000-0000	000-0000	000-0000	000-0000	000-0000	000-0000	000-0000	000-0000
Fine	000-0000	000-0000	000-0000	000-0000	-	-	-	-	-	-	-	-	-

	Flat End Cylinder	Long Pear	Round End Taper						Flat End Taper					
Shape No.	837	830	850	850	855	850	855	856	847	847	847	847	847	848
Size (1/10 mm)	012	014	012	014	016	018	018	018	012	014	016	018	018	023
Length/Tip (mm)	7.0	3.0	11/0.5	10/0.8	7/1.0	10/1.1	7/1.4	9/1.4	9/0.6	9/0.7	9/1.1	8/1.1	10/1.4	
Spring No.	515.7	365.4	770.11	771.10	772.7	772.10	773.7	773.9	700.9	701.9	702.9	703.8	703.10	
Supercoarse	000-0000	-	-	000-0000	000-0000	000-0000	000-0000	000-0000	-	-	000-0000	000-0000	000-0000	
Coarse	000-0000	000-0000	000-0000	000-0000	000-0000	-	000-0000	000-0000	000-0000	000-0000	000-0000	-	-	
Fine	-	-	000-0000	-	-	-	-	-	000-0000	-	-	-	-	

Turbo FG Diamonds (10 & 25 Packs Only)

	Turbo FG Diamonds														
Shape No.	368T	863T	836T	880T	880T	880T	850T	855T	850T	855T	856T	847T	847T	837KRT	850KRT
Size (1/10 mm)	023	014	012	010	012	014	014	016	018	018	018	016	018	012	018
Length (mm)	5.0	8.0	5.0	7.0	7.0	8.0	10.0	7.0	10.0	7.0	9.0	9.0	8.0	8.0	9.0
Spring No.	285.5T	265.8T	515.5T	KS0T	KS1T	KS2T	771.10T	772.7T	772.10T	773.7T	773.9T	702.9T	703.8T	837KRT012	850KRT018
Supercoarse	000-0000	000-0000	000-0000	000-0000	000-0000	000-0000	000-0000	000-0000	000-0000	000-0000	000-0000	000-0000	000-0000	000-0000	000-0000

	Turbo FG Diamonds					
Shape No.	878T	879T	862T	863T	862T	863T
Size (1/10 mm)	012	012	012	012	016	016
Length (mm)	8.0	10.0	8.0	10.0	8.0	10.0
Spring No.	250.8T	250.10T	260.8T	260.10T	267.8T	267.10T
Super Coarse	000-0000	000-0000	000-0000	000-0000	000-0000	000-0000

TO ORDER:
 For packs of 10 add
 -10 after order no.
 For packs of 25 add
 -25 after order no.

Regular FG Diamonds
 10 Per Package @ 00.00
 25 Per Package @ 00.00
Turbo Supercoarse FG Diamonds
 10 Per Package @ 00.00
 25 Per Package @ 00.00