



Power, Speed, Value. Spring Has It All!

High Quality Single-use Diamonds

Sold in Packs of 10 & 25

	Bulk Reduction Turbo		Occlusal Reduction & Finishing Wheel		Flame				Round		Flat End Taper					
Shape No.	837KRT	6055	368	909	879	862	863	863	801	801	847	847	847	847	848	
Size (1/10 mm)	012	018	023	035	018	012	012	014	014	023	012	014	016	018	023	
Length (mm)	8.0	8.0	5.0	1.1	9.0	8.0	10.0	8.0	1.4	2.3	9/0.6	9/0.7	9/1.1	8/1.1	10/1.4	
Spring No.	837KRT012	S18	285.5	860	256.9	260.8	260.10	265.8	120	135	700.9	701.9	702.9	703.8	703.10	
Super Coarse 10 Pack	-	-	138-9646	-	-	-	-	-	-	-	-	-	-	-	-	138-4845
Super Coarse 25 Pack	138-8856	138-4731	138-4597	138-1753	-	-	-	-	-	138-1990	-	-	138-8428	138-8578	138-8840	
Coarse 10 Pack	-	-	-	-	-	138-9281	138-6399	138-1048	-	-	-	-	-	-	-	
Coarse 25 Pack	-	-	-	-	138-7407	138-8864	138-0612	138-1103	138-1636	138-1770	138-5304	138-5553	138-5687	-	-	
Fine 10 Pack	-	-	-	-	-	-	138-4979	-	-	-	-	-	-	-	-	
Fine 25 Pack	-	-	138-4323	-	-	138-8510	138-8995	-	-	-	138-5033	-	-	-	-	

Sold in Packs of 10 & 25

	Beveled Cylinder	Flat End Cylinder			Round End Taper					Turbo	
Shape No.	878	835	836	837	850	850	855	850	855	856	856T
Size (1/10 mm)	012	008	012	012	012	014	016	018	018	018	018
Length (mm)	8.0	4.0	5.0	7.0	11/0.5	10/0.8	7/1.0	10/1.1	7/1.4	9/1.4	9/1.4
Spring No.	250.8	514.4	515.5	515.7	770.11	771.10	772.7	772.10	773.7	773.9	773.9T
Super Coarse 10 Pack	-	-	-	-	-	-	138-1574	138-9851	-	138-9424	-
Super Coarse 25 Pack	138-6779	-	138-4186	138-4521	-	-	138-0160	138-0396	138-0795	138-1276	138-2381
Coarse 10 Pack	-	138-3823	-	-	-	138-6939	138-2881	-	-	138-2991	-
Coarse 25 Pack	138-6655	-	138-3954	138-4288	138-9362	138-9764	138-0011	-	138-0661	138-1018	-

PRICING

TURBO DIAMONDS
25 Per Package @ 72.49
= **2.92** per diamond.

REGULAR DIAMONDS
10 Per Package @ 15.99 =
1.60 per diamond.
25 Per Package @ 35.29 =
1.41 per diamond.