

# NEW!!! Spring L.E.D. Headlamp



**\*1-YEAR WARRANTY**

## Features

- *Brightest in the Industry*
- *Dual lights to increase viewing and reduce shadowing*
- *Quick Batt. Replacement, 2 batteries included.*
- *Battery pack can be placed on headset or belt clip.*
- *Headset is fully adjustable for comfort*
- *Simple durable design*



Available through your major dental dealer.

## Specifications

- *3-4 hour battery life*
- *Recharge time 3-4 hours per battery*
- *Weight: 10 oz.*
- *Weight without battery (battery located on belt) 6 oz.*
- *Bulb Life 5000hrs.*
- *27,000 LUX/ 1 ft.*
- *6500K correlated color temp.*



PATENT PENDING