



SPRING HIGH SPEED HANDPIECE FOR KAVO COUPLER



Key Features:

- Ceramic Bearings (Made in the USA): Lighter, Smoother Therefore Less Friction, Wear and Less Heat
- Extremely Low Noise Level
- Patented Push Button Chuck for Reliability and Security
- Evenly Balanced for Comfort and Ergonomic Shape
- Fiber Optic Push Button Swivel
- Glass Rod for Optics – Does Not Yellow
- Length: 4.5 In = 11.4 Cm
- Made in the USA
- Shell One- Year Warranty: Turbine Six-Month Warranty
- Speed 400,000 Rpm
- Torque 18 Watts
- Precise Triple Port Water Spray Avoids Overheating

Options:

- LED Light Available for Kavo Type Coupler