


the CURE



- High Output – 6 seconds at high intensity for most composites
- Shaped to be held like a handpiece
- 470 nanometers
- Low and High Intensity (500-1000mW/cm²)
- Size: 6 inches
- Weight: 1 ounce

**\$499
retail**

At first glance you might think that this is nothing but an elegantly designed light. You would be right of course, but that is just the beginning!

FACT #1

No other halogen or LED has higher power, and the peak power is at the wavelength of 470 nanometers, perfect for composites. So curing time is faster than all other lights.

FACT #2

It's better than a cordless. The charger for a cordless occupies counter space and has the power cord draped from the counter to the wall socket.

the CURE will never need a counter top for storage. With holders supplied it can be mounted anywhere on any unit, wall or cabinet. Your choice!

FACT #3

the CURE is the only light with two on/off buttons. They are opposite each other, so that you can cure uppers or lowers without a tortured finger position.

FACT #4

the CURE is the only light which allows you the option of step curing.

You can switch from low to high power as you are curing.

FACT #5

If you think those facts are impressive wait until you use it, it's so light its like lifting air.

FACT #6

The price? It could be the highest. It is the most advanced. But, the price is only \$499.00. We keep it low so that you'll come back for your second and third operator. And you can try it with the comfort of knowing that for 30 days you can return it for full money back!

